



Ferrous Iron Reagent Powder Pillows, pk/100

Product #: 103769
USD Price: Contact Hach
Available

For determination of ferrous iron by the 1,10-phenanthroline method. Hach method 8146. Range: 0.02-3.00 mg/L Fe. Pack of 100 powder pillows.

Specifications

Instrument:	All except PC II
Method:	8146
Method Name:	1, 10 Phenanthroline
Number of tests:	100
Parameter:	Iron, Ferrous
Platform :	Powder Pillows
Quantity:	100 pk
Range:	0.02 - 3.00 mg/L Fe ²⁺

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)