



Manganous Sulfate Powder Pillows, pk/50

Product #:	107166
USD Price:	Contact Hach
Available	

For determination of Dissolved Oxygen by the Winkler titration method. Range: 1 to >10 mg/L. Hach methods 8215 and 8229. Sample size: 300 mL.

Specifications

Number of tests:	50
Parameter:	Dissolved Oxygen
Platform :	Powder Pillows
Range:	1 - greater than 10 mg/L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)