



Triazole Reagent Powder Pillows, pk/100

Product #:	2141299
USD Price:	Contact Hach
Available	

For the determination of benzotriazole (1-16 mg/L) or tolyltriazole (1-20 mg/L) by the UV Photolysis method. Hach method 8079. Requires the use of a UV Lamp (Cat. No. 2082800 or 2082802). Pack of 100 powder pillows

Specifications

Instrument:	All except PC II, DR 900
Method:	8079
Method Name:	UV Photolysis
Number of tests:	100
Parameter:	Benzotriazole, Tolyltriazole
Platform :	Powder Pillows
Quantity:	pk
Range:	1.0 - 20.0 mg/L Tolyltriazole
Range 2:	1.0 - 16.0 mg/L Benzotriazole

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)