

Ultra-low temperature refrigerated vapor trap (-105°C)

Condense and trap a wide variety of solvent vapors that can contaminate and deteriorate vacuum pump oil with Thermo Scientific™ Savant Refrigerated Steam™ Traps.

- Rapid cooling minimizes waiting time
- Easy-to-remove glass bottle for quick solvent changes
- Fast cooling to operating temperature
- Easy-to-use glass condensing bottle allows liquid to be collected in a separate jug for easy draining
- Hermetic cooling system saves consumable refrigerants such as dry ice and acetone

