

FLN456LED-AVL

Ensure optimal safety for critical samples by using Thermo Scientific™ CO_2 and LN_2 Support Systems for ULT Freezers, for use with Thermo Scientific™ TSX, Forma™, HERAfreeze™, Revco™, and TS Series freezers.

Features:

- LN_2 standby system maintains liquid nitrogen temperature at -80°C with liquid nitrogen
- CO_2 reserve system that maintains temperature at -67°C with CO_2

Ordering Information:

Available in factory-installed or site-installed format. Models for a 300 box capacity are independent. The installation of on-site units must be carried out by qualified personnel

