

Adpt 250mL to 10mL Pyrex EA

- Combine the Thermo Scientific™ Fiberlite™ F14-6 x 250LE rotor adapters with the Fiberlite F14-6 x 250LE fixed-angle rotor.
- Attention:
- Adapters 00369 and 00456 should be used together for 7 mL round-bottom tubes
- Ordering Information:
- The Fiberlite F14-6x250LE (75003662) fixed-angle rotor body is also available. Contact your sales representative for more information.

Tube Size 50 mL
Tube Type Tapered
Capacity 50 mL
Type Rotor adapter
Unit Size Pack of 2

