

Free standing field installed LN2 backup system 230V/50Hz with British plug fits UxF300V HFU300TV 88300V TSU300V

- Ensure optimal safety for critical samples by using Thermo Scientific™_{CO2} and_{LN2} Support Systems for ULT Freezers, for use with Thermo Scientific™ TSX, Forma™, HERAfreeze™, Revco™, and TS Series freezers.
- Features:
 - LN₂ standby system maintains liquid nitrogen temperature at -80°C with liquid nitrogen
 - CO₂ reserve system that maintains temperature at -67°C with_{CO2}
- Ordering Information:
 - Available in factory-installed or site-installed format. Models for a 300 box capacity are independent. The installation of on-site units must be carried out by qualified personnel.

