

**Free Standing 115V/60Hz field installed CO2 backup
system;115V/60Hz; for use with RLE30086A; HLE30086A;
8930086A; TLE30086A**

- Ensure optimal safety for critical samples by using Thermo Scientific™ CO₂ and LN₂ Support Systems for ULT Freezers, for use with Thermo Scientific™ TSX, Forma™, HERAfreeze™, Revco™, and TS Series freezers.
- Features:
 - LN₂ standby system maintains liquid nitrogen temperature at -80°C with liquid nitrogen
 - CO₂ reserve system that maintains temperature at -67°C with CO₂
- Ordering Information:
 - Available in factory-installed or site-installed format. Models for a 300 box capacity are independent. The installation of on-site units must be carried out by qualified personnel.

