

Free Standing 115V/60Hz field installed LN2 backup system 115V/60Hz for use with UxF300A HFU300TA 88300A TSU300A

- Ensure optimal safety for critical samples by using Thermo Scientific™ CO₂ and LN₂ Support Systems for ULT Freezers, for use with Thermo Scientific™ TSX, Forma™, HERAfreeze™, Revco™, and TS Series freezers.
- Features:
 - LN₂ standby system maintains liquid nitrogen temperature at -80°C with liquid nitrogen
 - CO₂ reserve system that maintains temperature at -67°C with CO₂
- Ordering Information:
 - Available in factory-installed or site-installed format. Models for a 300 box capacity are independent. The installation of on-site units must be carried out by qualified personnel.

