



Low-volume flow chamber, PVDF, 3422-series sensors

Product #: 1000G3316-101

USD Price: Contact Hach

Available

For use with the 3422 Contacting Conductivity Sensors K=0.05 Cell Constant Only, Kynar (PVDF) 1/2-Inch Pipe Tee, Limits sample volume to approximately 20 ml for high-purity water applications.

Specifications

Material: Kynar® (PVDF)

Process Connection: 0.5 "

Special notes: These tees limit

sample volume to approximately 20 mL for high-purity water applications.

Special Notes: For use only with a 0.05 cell constant sensor.