



Molybdate Reagent Powder Pillows for HR Silica, 25mL, pk100

Product #:	104199
USD Price:	Contact Hach
Available	

For determination of high range silica by the Silicomolybdate method. Hach method 8185. Sample size: 25 mL. Pack of 100 powder pillows

Specifications

Method:	8185
Method Name:	Silicomolybdate
Number of tests:	100
Parameter:	Silica
Platform :	Powder Pillows
Quantity:	100 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)