



## Acid Reagent Powder Pillows, 25 mL, pk/100

Product #:

104299

USD Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of high range silica by the silicomolybdate method (Hach method 8185) Range: 1-100.0 mg/L SiO<sub>2</sub>. Also used in the iodometric method (Hach method 10100) for the determination of high range hypchlorite (50-150 g/L as Cl<sub>2</sub>). Pack of 100 powder pillows.

---

### Specifications

|                  |                 |
|------------------|-----------------|
| Method:          | 8185            |
| Method Name:     | Silicomolybdate |
| Number of tests: | 100             |
| Parameter:       | Silica          |
| Platform :       | Powder Pillows  |
| Quantity:        | 100 /pk         |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)