



Beaker, polypropylene, low form, 400 mL, 6/pk

Ref. producto: 108078
EUR Precio: Contáctenos
Envío dentro de las próximas 2 semanas

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Especificaciones

Capacidad:	400 mL
Diámetro:	99.1 mm with spout
Material carcasa:	Polypropylene
Volumen / tamaño de embalaje:	6 /pk