



Cylinder, Graduated, Polypropylene, 250 mL $\pm$ 1.4 mL, 2.0 mL divisions

Product #:	108146
USD Price:	Contact Hach
Available	

Polypropylene construction of both cylinder and base. All sizes have generous pour spout, large molded-in, and easy-to-read graduations. Can be chemically sterilized without affecting accuracy. Autoclavable, but accuracy will be adversely affected. Cylinders meet accuracy requirements of ASTM Class B, E1272, ISO Standard 6706, and CFR 21, Part 177.1520

Specifications

Body material:	PP
Capacity:	250 mL
Certifications:	None
I.D.:	34.5 mm
Length:	353 mm
Measurement Method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric