



Cylinder, Graduated, Polypropylene, 500 mL $\pm$ 2.6 mL, 5.0 mL divisions

Product #:	108149
USD Price:	Contact Hach
Available	

Polypropylene construction of both cylinder and base. All sizes have generous pour spout, large molded-in, and easy-to-read graduations. Can be chemically sterilized without affecting accuracy. Autoclavable, but accuracy will be adversely affected. Cylinders meet accuracy requirements of ASTM Class B, E1272, ISO Standard 6706, and CFR 21, Part 177.1520

Specifications

Body material:	PP
Capacity:	500 mL
Certifications:	None
Instrument:	All except DR/4000, PC II
Measurement Method:	To Deliver/To Contain
Method:	8006
Quantity:	1 each
Type:	Single Metric