



Flask, volumetric, polypropylene, 250 mL

Ref. producto: 1406046
EUR Precio: Contáctenos
Llamar para consultar fecha de envío

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Especificaciones

Capacidad:	250 mL
Volumen / tamaño de embalaje:	1 each