



## Flask, Volumetric, Polypropylene, 1000 mL

|            |              |
|------------|--------------|
| Product #: | 1406053      |
| USD Price: | Contact Hach |
| Available  |              |

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

---

### Specifications

|           |         |
|-----------|---------|
| Capacity: | 1000 mL |
| Quantity: | 1 each  |