



Iron Color Disc Test Kit, Model IR-10

Product #: 146430
USD Price: Contact Hach
Call for ship date

Using the simple color disc method, this kit measures Iron (medium range: 0.1 - 10 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Quickly and accurately measure iron in the field or in the lab

Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Method Name:	Color disc/ 1,10 Phenanthroline
Model:	IR-10
Number of tests:	100
Parameter:	Iron, medium range - as Fe
Platform :	Color Disc
Range:	0 - 10 mg/L
Smallest Increment:	0.1

What's in the box?

Includes one reagent, two color discs, comparator box, two viewing tubes, instructions and carrying case.