



Balance, Hanging Pan, 300 g

Product #: 1476000
USD Price: Contact Hach
Call for ship date

Ohaus Cent-O-Gram: good, general purpose hanging pan balance with 300 gram weighing capacity, sensitive to 0.01 g. Readings are taken from four beams of balance. Smallest increment is 0.01 g. Stainless steel weighing pan with handle and two pouring spouts.

Pan diameter: 9.25 cm (3.64"). Depth: 1.3 cm (0.5"). Shipping weight: 2.1 kg (4.8 lb).

Economical

Perfect for general laboratory weighing.

Accurate

Unit has magnetic-damping, self-storing specific gravity stand, and zero adjustment knob.

Specifications

Accuracy 2:	± 0.01 g Readability
Capacity:	311 g
Communication Capabilities:	None
Power:	None
Quantity:	1 each