



## Free Amino Nitrogen Reagent Set

Product #:

45950

USD Price:

Contact Hach

Call for ship date



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

### Simplify your fermentation efficiency by monitoring your Free Amino Nitrogen.

For use with your Hach DR6000 spectrophotometer and designed to simplify brewing analysis.

Testing for FAN is a part of standard brewing quality control analyses. This enables the protein content to be estimated and serves as an indicator of beer quality.

---

### Specifications

Instrument:

DR6000

Parameter:

Free Amino Nitrogen

---

### What's in the box?

Di-Sodium-Hydrogen Phosphate, 500 g (170234) Potassium Phosphate, Monobasic, ACS, 454 g (17001H) Ninhydrin, 99%, 10 g (129523) D-Fructose, 99%, 250 g (135129) Potassium Iodate, ACS, 100 g (606932) Glycine, ACS, 100 g (2111132)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)