



Tensette® Pipet, 0.1-1.0 mL

Product #:	1970001
USD Price:	Contact Hach
Available	

Perform routine laboratory pipetting quickly, economically, and accurately with Hach's automatic TenSette® Pipet. Two models dispense from 0.1 to 1 ml or from 1 to 10 ml in 0.1 or 1.0 ml increments, respectively. The pipet saves time, reduces glassware cleanup and assures accuracy. The 0.1 to 1.0 mL model is specifically designed for use with Hach Voluette and PourRite standards in Standard Additions accuracy checks. Set includes pipet, manual and 100 disposable tips.

Specifications

Quantity:	1 each
-----------	--------