



Carboy, Lowboy, 15-L Capacity

Product #:	2075460
USD Price:	Contact Hach
Available	

Made of high-density polyethylene. Low-profile design fits into narrow spaces, on shelves, and in refrigerators. Stacks up to three high to save space. Large, integral, hand grip makes carrying easy. 8 L size can be stacked on 15 L size. For storage and dispensing of reagents, distilled water, titrants, and other liquid storage needs. Will tolerate strong acids and bases, alcohols, and ketones. Not recommended for use with other organic solvents. Leak-proof. Includes spigot.

Specifications

Capacity:	15 L
Depth:	14.4 Inches
Dimensions: Depth:	36.5 cm
Length:	201 mm
Quantity:	1 each
Width 2:	14.5 Inches