



Demineralizer Cartridge, UltraPure

Product #: 2079100
USD Price: Contact Hach
Ships within 2 weeks

Color change indicator identifies when resin is exhausted. Largely free of silica and carbon dioxide. Designed for non-pressure applications, cartridge includes 3/8 inch barb hose connections on each end.

Ultra Pure:

No dye in resin. Largely free of silica and carbon dioxide.

Organic Removal:

Activated carbon removes organic and chlorines from feed water. Place ahead of Ultra Pure or High Capacity cartridge. Replace after second demineralizer cartridge replacement.

High Capacity:

Lower purity than Ultra Pure.

Color Indicator:

Ultrapure mixed bed ion exchange changes from dark purple to a carmel color as it exhausts.

Specifications

Quantity: 1 each