



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 500 mL

Ref. producto: 2089679
EUR Precio: Contáctenos
Llamar para consultar fecha de envío

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Especificaciones

Capacidad: 500 mL
Volumen / tamaño de embalaje: 6 /pk