



# Pipet, Serological, Sterile, Polystyrene, 2 mL, White, 1000/pk

|                      |              |
|----------------------|--------------|
| Product #:           | 2092656      |
| USD Price:           | Contact Hach |
| Ships within 2 weeks |              |

Sterile, disposable plastic pipet. Individually wrapped. Color coded. Calibrated “To Deliver”. Complies with ASTM E934.

## Specifications

|                       |          |
|-----------------------|----------|
| Capacity:             | 2 mL     |
| Negative Graduations: | 0.3 mL   |
| Quantity:             | 1000 /pk |
| Resolution:           | 0.01 mL  |