



Grab sampler

Ref. producto: 2095200
EUR Precio: Contáctenos
Llamar para consultar fecha de envío

Guarantees underwater sampling without surface contamination. Assists in conforming to procedures described in the USEPA Handbook for Monitoring Industrial Wastewater. Clamp the sterilized, capped bottle onto the grab sampler before sterilizing. Remove cap, collect sample, and replace cap without touching the bottle or getting wet. Sampler is constructed of 1.9 cm (0.75") rigid square aluminum tubing with baked enamel finish. Comes complete with clamp and graduated 1000 mL borosilicate glass sampling bottle with PTFE-lined cap. Adjustable clamp for smaller bottles also provided. Sampler length: 1.8 m (6'). Weight: 1.8 kg (4 lb).

Especificaciones

Volumen / tamaño de embalaje: 1 each