



**Be Right™**



## Brush, Camel's Hair

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2103200</b> |
| USD Price:        | Contact Hach   |
| Available         |                |

For cleaning delicate instruments or to ensure complete transfer of samples when weighing. Soft, brown hair is cemented into a metal holder. Hairs are cut to the same length—14 mm (0.5"). Can be attached to Micro Spatula, Prod. No. 1225600.

---

### Specifications

|           |        |
|-----------|--------|
| Quantity: | 1 each |
|-----------|--------|