



# Pipet, Mohr, Polypropylene, 5 mL, 12/pk

|                      |              |
|----------------------|--------------|
| Product #:           | 210667       |
| USD Price:           | Contact Hach |
| Ships within 2 weeks |              |

Graduated, translucent polypropylene pipets. Ideal for silica analysis. Autoclavable (Temperatures above 60°C can cause changes to volumes).

## Specifications

|             |        |
|-------------|--------|
| Capacity:   | 5 mL   |
| Quantity:   | 12 /pk |
| Resolution: | 0.1 mL |