



Cylinder, graduated, polymethylpentene, 250 mL $\pm$ 1.4 mL, 2.0 mL divisions

Ref. producto: 217246
EUR Precio: Contáctenos
Envío dentro de las próximas 2 semanas

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated "To Contain/To Deliver" at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Especificaciones

Capacidad:	250 mL
Certificaciones:	Class B
Material carcasa:	PMP
Método de medición:	To Deliver/To Contain
Tipo:	Métrica simple
Volumen / tamaño de embalaje:	1 each