



Cylinder, Graduated, Polymethylpentene, 500 mL $\pm$ 2.6 mL, 5.0 mL divisions

Product #:	217249
USD Price:	Contact Hach
Available	

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated “To Contain/To Deliver” at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Specifications

Body material:	PMP
Capacity:	500 mL
Certifications:	Class B
Measurement Method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric