



Lamp, Ultraviolet, Long-Wave 115 Vac, 60 Hz

Product #:

2184300

Contains mercury. Dispose of according to local, state and federal regulations.

USD Price:
Available

Contact Hach

Mid-size (6 Watt) high intensity, hand-held and line-operated. Intensity at 15.2 cm (6"), 750 $\mu\text{W}/\text{cm}^2$. Filters: Two 6 x 7.6 cm (2.4 x 3") in snap on/off frame. Overall dimensions: 37.5 x 8.9 x 6.7 cm (14.8 x 3.5 x 2.6"). Weight: 0.91 kg (2 lbs). UL listed. Meets existing OSHA electrical standards. 365 nm.

The Fastest Presence/Absence Test on the Market at 16 Hours

Fastest EPA-approved test for coliforms and E. coli. Test as

late as 5:00 PM and still have your results by 9:00 AM the next day.

One Product that Supports 16- or 22-hour Test Methods

Allows for a "standard" test that requires less set-up, or a

"rapid" test when necessary. One media for both tests reduces

complexity in certification and inventory.

Results are Valid for 32 Hours

Saturday's tests can be read on Monday—no need to work all

weekend. Results stay valid for 800% longer than the competition.

Specifications

Quantity: 1