



Distillation Apparatus Set

Product #:	2265300
USD Price:	Contact Hach
Call for ship date	

Fast set-up for sample pretreatment.

Fast set-up for sample pretreatment. This Distillation set includes EPA approved glassware ONLY. Easy to set up because you only receive the parts you need. Illustrated, step-by-step procedures simplify distillation.

***Heater and support not included.**

Use in EPA distillations

Distillation sets include everything you need to satisfy EPA distillation requirements

Complete glassware sets available

Easy to set up because you only receive the parts you need

Easy to follow instructions

Illustrated, step-by-step procedures simplify distillation

Versatile

Arsenic and Cyanide add-on accessories available

Specifications

Quantity:	Each 1
-----------	--------

What's in the box?

Includes distillation flask, J-tube, condenser, drip tube, erlenmeyerflask, and non-Hg thermometer.