



Beaker, Berzelius glass, Tall Form, without spout, 200mL 12/pk

Product #: 2276175
USD Price: Contact Hach
Ships within 1 week

200 mL Berzelius style (tall form) beakers from Kimax, are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Glass
Capacity:	200 mL
Package size UOM:	per pack
Quantity:	12 /pk