



Solución tampón de pH, pH 10,01 (500 mL)

de producto: 2283649
COP Precio: Contáctenos
Disponible

pH Buffer standard solutions for laboratory and field analyses

Solución tampón de calidad, codificada en color azul, certificada por el NIST.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Especificaciones

Ajustes del envase:	Bottle
Concentración:	10,01 pH $\pm$ 0,02 pH a 25 °C
Condiciones de almacenamiento:	10 - 25 °C (protect from light)
Descripción:	Solución tampón de pH, codificada por colores, pH 10,01 (azul)
Parámetro:	pH
Trazable a NIST:	Sí
Valor(es) de pH:	10,01
Vida útil:	12 meses
Volumen / tamaño de embalaje:	500 mL
¿Codificada por colores?:	Sí