



sensION Sodium Ion Selective Electrode (ISE) Analysis Package

Product #:

2348100

USD Price:

Contact Hach

Call for ship date

sensION Sodium Analysis Package includes Sodium Platinum Series Combination Electrode, two ammonium chloride reference cartridges, 100 and 1000 mg/L standards, reagents, beaker, and manual. The sensION Sodium ISE is compatible with Hach sensION, sensION+, and H-Series systems.

Everything you need for ISE measurements in one package.

Hach's ISE Analysis Packages provide a complete solution for ion selective measurements, including ISE, standard solutions, ionic strength adjustor, and replacement electrolyte cartridges.

Maximum life through reduced maintenance.

sensION Platinum Series probes' replaceable gel reference electrolyte cartridge is a fast and easy maintenance solution that maximizes the life of your electrode.

Accurate

You don't just need an answer, but the right one, and fast.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Specifications

Cable Length:	1 m
Diameter:	12 mm
Filling Solution:	2595802
ISA Required:	4451569
Junction:	Open
Length:	240 mm
Material:	Sensor Body: TPX
Probe Type:	Laboratory combination
Range:	0.01 - saturation mg/L Na ⁺
Sensor Type:	Glass
Special Feature:	Platinum Series
Temperature Range:	0 - 50 °C
Temperature Range 2:	N/A N/A N/A

Thermistor: -

What's in the box?

Kit includes: sensION 51925 combination Platinum Series Sodium ISE (no. 5192500), two replacement electrolyte cartridges (no. 2595802), two Hach sodium standard solutions (100 mg/L (1 L, no. 2318153) and 1000 mg/L (500 mL, no. 1474949)), Sodium Chloride (454 gr., no. 18201H), Hach sodium ionic strength adjustor (ISA) Powder Pillows (pack of 100, no. 4451569), 50-mL beaker (no. 108041), and ISE User Manual.

Required Accessories

- Ammonium Chloride Reference Cartridge, pk/2 (Item 2595802)
- Sodium/Potassium Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 4451569)