

ACIDO LACTICO 85% RA ACS

About this item

Boiling Pt:122 °C (20 hPa)

Density:1.20 g/cm³

Storage Temperature:Ambient

UN:3265

ADR:8,II :90.08 g/mol

EINECS:200-018-0 :GHS05

CAS Number:50-21-5 :598-82-3

Formula:H₃CCH(OH)COOH

MDL Number:MFCD00004520

Specifications

Pack type:	Glass bottle, amber, narrow mouth
Pk:	500 mL

Specification Test Results

Meets ACS Reagent Chemical Requirements:	85.0 - 90.0%
For Laboratory, Research, or Manufacturing Use:	≤0.02%
ACS - Assay (as CH ₃ CHOHCOOH):	85.0 - 90.0 %
ACS - Residue after Ignition:	≤ 0.02 %

Seq:	2725
------	------

ACS - Chloride (Cl):	≤ 0.001 %
ACS - Sulfate (SO ₄):	≤ 0.002 %
ACS - Substances Darkened by H ₂ SO ₄ :	Passes Test
Trace Impurities - ACS - Heavy Metals (as Pb):	≤ 5 ppm
Trace Impurities - Iron (Fe):	≤ 5 ppm