

About this item

CAS Number:10043-52-4

Formula:CaCl₂

MDL Number:MFCD00010903

Boiling Pt:1600 °C (1013 hPa)

Melting Pt:772 °C

Density:2.174 g/cm³ (20 °C)

:Management of Change (MOC) category = N

:110.98 g/mol

Storage Temperature:Ambient

EINECS:233-140-8

REACH:01-2119494219-28

:GHS0

Specification Test Results

Assay (CaCl ₂) (by Ag titrn):	≥96.0 %
Titrate Base (meq/g):	≤0.006

Mesh - On U.S. No. 20 Sieve:	≤10 %
Mesh - Thru U.S. No. 200 Sieve:	≤30 %