

About this item

CAS Number:1305-78-8

MDL Number:MFCD00010911:56.08 g/mol

Boiling Pt:2850 °C (1013 hPa)

Density:3,38 g/cm³ (20 °C)

:GHS05

Formula:CaO

Melting Pt:2580 °C

EINECS:215-138-9

Storage Temperature:Ambient

Specifications

Pack type:	Poly bottle, blue T/E cap
Pk:	500 g

Seq:	3645
------	------

Specification Test Results

For Laboratory, Research, or Manufacturing Use:	≥98.0%
Assay (BaCl ₂):	≥98.0 %

Insoluble Matter:	≤0.02 %
pH of 5% Solution at 25°C:	4.5 - 7.5
Heavy Metals (as Pb):	≤0.002 %