

About this item

Melting Pt:300 °C

Storage Temperature:Ambient

Formula: $\text{H}_3\text{C}(\text{CH}_2)_5\text{CH}_2\text{SO}_3\text{Na}$

EINECS:245-210-5 :GHS07

CAS Number:22767-50-6

MDL Number:MFCD00007543

Specifications

Pack type:	Glass bottle, amber, wide mouth
Pk:	25 g

Seq:	3308
------	------

Specification Test Results

For Laboratory, Research, or Manufacturing Use:	≥98.0%
Assay (as $\text{C}_7\text{H}_{15}\text{SO}_3\text{Na}$) (by Acidimetry):	≥98.0%
UV Absorbance of 0.25 M Solution - 250 nm:	≤0.050
UV Absorbance of 0.25 M Solution - 240 nm:	≤0.050

UV Absorbance of 0.25 M Solution - 230 nm:	≤0.050
UV Absorbance of 0.25 M Solution - 220 nm:	≤0.060
UV Absorbance of 0.25 M Solution - 210 nm:	≤0.080
UV Absorbance of 0.25 M Solution - 200 nm:	≤0.200