

PERÓXIDO HIDRÓGENO 30%

About this item

Boiling Pt:~107 °C (1013 hPa) :GHS03

Melting Pt:< 0 °C :GHS05

Density:1.11 g/cm³ (20 °C) :GHS06

UN:2014 :GHS08

ADR:5.1 (8),II:Management of Change (MOC) category = N :34.01

EINECS:231-765-0

Storage Temperature:Ambient

CAS Number:7722-84-1

Formula:H₂O₂

MDL Number:MFCD00011333

Specifications

Pack type:	Poly bottle
Size:	500 mL

Cat. No.:	JT2189-1
Supplier no.:	2189-01

Specification Test Results

Assay (H ₂ O ₂):	29.0 - 32.0 %
Residue after Evaporation:	≤0.0300 %
Nitrogen Compounds (as N):	≤0.001 %
Phosphate (PO ₄):	≤0.02 %

Sulfate (SO ₄):	≤0.001 %
Chloride (Cl):	≤5 ppm
ACS - Heavy Metals (as Pb):	≤1 ppm
Trace Impurities - Iron (Fe):	≤0.5 ppm