

Biftalato de Potasio Crist.

About this item

CAS Number:877-24-7

Formula:HOOC₆H₄COOK

MDL Number:MFCD00013070 :204.22 g/mol

Melting Pt:295...300 °C

Density:1.6362 g/cm³ (20 °C)

:Management of Change (MOC) category = N

Storage Temperature:Ambient

EINECS:212-889-4

Specification Test Results

Assay (1-KOCC ₆ H ₄ -2-COOH) (dried basis):	99.95 - 100.05 %
Insoluble Matter:	≤0.005 %
pH of 0.05M Solution at 25°C:	4.00 - 4.02
Chlorinated Compounds (as Cl):	≤0.003 %

Sulfur Compounds (as S):	≤0.002 %
Sodium (Na):	≤0.005 %
Trace Impurities - ACS - Heavy Metals (as Pb):	≤5 ppm
Trace Impurities - Iron (Fe):	≤5 ppm