

Bromato de Potasio RA.

About this item

Melting Pt:350 °C

Density:3.218 g/cm³ (20 °C)

Storage Temperature:Ambient

EINECS:231-829-8

UN:1484

ADR:5.1,II

:Management of Change (MOC) category = N

CAS Number:7758-01-2

Formula:KBrO₃

MDL Number:MFCD00011359

Specifications

Pack type:	Glass bottle
Size:	500 g

Cat. No.:	JT2992-1
Supplier no.:	2992-01

Specification Test Results

Assay (KBrO₃) (by iodometry):	≥99.8 %
Insoluble Matter:	≤0.005 %
pH of 5% Solution at 25°C:	5.0 - 9.0
Bromide (Br):	Passes Test

Sulfate (SO₄):	≤0.005 %
Iron (Fe):	≤0.002 %
Sodium (Na):	≤0.010 %
Trace Impurities - ACS - Heavy Metals (as Pb):	≤5 ppm

