

Carbonato de Sodio Monohid.

About this item

CAS Number: 5968-11-6

Formula: $\text{Na}_2\text{CO}_3 \cdot \text{H}_2\text{O}$

MDL Number: MFCD00149177

Melting Pt: 107.85 °C

Molecular Weight: 124.00 g/mol

Density: 2.25 g/cm³ (20 °C)

: Management of Change (MOC)
category = R

Storage Temperature: Ambient

EINECS: 207-838-8

Specification Test Results

NF - Identification A:	Passes Test
NF - Identification B:	Passes Test
NF - Water (H ₂ O):	12.0 - 15.0 %
NF - Assay (Na ₂ CO ₃) (anhydrous basis):	99.5 - 100.5 %
FCC - Identification A:	Passes Test

FCC - Identification B:	Passes Test
FCC - Lead (Pb):	≤4 mg/kg
FCC - Loss on Drying:	12.0 - 15.0 %
FCC - Assay (Na ₂ CO ₃) (anhydrous basis):	99.5 - 100.5 %