

Tiras indicadoras del pH pH 2,5 - 4,5

Descripción	
Número de catálogo	109541
Replaces	9581-1; 9581-3; 9581
Visión de conjunto	The MQuant® pH tests are suited for all media in environmental analysis and in industrial in-process controls. Our special method of manufacture gives the pH indicator strips decisive advantages over the indicator paper: since the indicator does not bleed, the strip can be left in the measurement medium without contaminating it. Using these test strips is fairly easy, the pH value is measured semiquantitatively by visual comparison of the reaction zone on the test strip with the fields of the color scale, perfectly matching the reaction color. The clever Safetyedge boxes allow for safe storage and convenient withdrawal of the test strips. For tests for other parameters that require a specific pH range of the sample (e.g. Spectroquant®, see instructions of use), we recommend the use of MQuant® pH test strips as pre-test.
Información del producto	
Código HS	3822 00 00
Quality Level	MQ100
Información de seguridad según el GHS	
Clase de almacenamiento	10 - 13 Otros líquidos y sustancias sólidas
WGK	NWG ningún peligro para el agua

Información de almacenamiento y transporte

Almacenamiento	Almacenar entre +15°C y +25°C.
-----------------------	--------------------------------