

Módulo de pretratamiento Progard S2

Descripción	
Número de catálogo	PR0G0T0S2
Descripción	Módulo de pretratamiento Progard S2
Visión de conjunto	The Progard® Pretreatment Pack has been designed to pretrat water by removing large particles, organic compounds and free chlorine from tap water. This protects the reverse osmosis (RO) membrane against mineral scaling, organic fouling and chlorine oxidation. Benefits: • Optimal protection of reverse osmosis (RO) membrane • Ensures system lifetime and performance are optimized
Información del producto	
Componentes	Polyphosphate for anti-scaling protection Silver-impregnated activated carbon for free chlorine reduction and bacteriostatic protection Depth-filtration for particles removal e-Sure tag RFID technology for safety and automatic pack parameters archiving
Aplicaciones	
Aplicación	For optimal protection of Reverse Osmosis membranes
Dimensiones	
Altura	39.5 cm (15.5 in)
Información sobre embalaje	
Cantidad	1