



Legion BLA3

Non-sterile Disposable Latex
Cleanroom Glove

Non-sterile latex cleanroom gloves, ensuring elbow length-protection and easy double-donning

- **Heightened flexibility and comfort:** BioClean™ Legion BLA3 cleanroom gloves' soft, stretchy formulation offer a perfect blend of flexibility and comfort
- **Added protection:** These ambidextrous cleanroom gloves are designed to provide coverage up to the elbow, granting wearers hand, forearm and wrist protection and greatly reducing cross contamination risks
- **Reduced contamination risks:** BioClean™ Extra BLA3 non-sterile gloves come in non-particulating EasyTear packaging, for further reduced contamination risks
- **Enhanced strength:** These powder-free gloves' textured surface enhances grip performance, while a beaded cuff bolsters strength
- **Increased practicality:** These BioClean™ gloves' robust formulation also makes them ideal for quick, trouble-free double-donning, while their ambidextrous nature ensures added convenience



KEY FEATURES & BENEFITS

- **Soft, stretchy film:** For flexible, comfortable cleanroom gloves
- **Elbow-length design:** Greater coverage, less cross contamination
- **Non-particulating EasyTear packaging:** Lesser contamination risks

Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D

WWW.ANSELL.COM

Ansell



Legion BLA3

Non-sterile Disposable Latex
Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Natural Rubber Latex
Color	Natural
Shape	Ambidextrous
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory	CE 0493, EN 421:2010, EN 455 Part 2, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, Category III, UKCA
Packaging	100 pieces per sealed inner PE bag; one inner PE bag per sealed outer PE bag; 10 outer bags per lined carton (1000 pieces). More sustainable packaging: These gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Malaysia
Available sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Chlorinated
Cleanroom Class	Class 100/ISO 5
Shelf Life	Five (5) years from date of manufacture.
Tested for use with Chemotherapy Drugs	No
Protein Level	50 µg/g or less of total extractable protein
Anti-static	No



Legion BLA3

Non-sterile Disposable Latex
Cleanroom Glove

Physical Properties					Testing Method
Sizes	S	M	L	XL	EN ISO 21420
Typical Length (mm/in)	400 / 16				
Palm Width (mm/in)	83/3.3	95/3.7	106/4.2	115/4.5	
Freedom from Holes	1.5 AQL Performance				EN 374-2
Typical Particle Count ≥0.5µm (counts / cm²)	<1500				IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.17 / 6.69				EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.20 / 7.87				EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.11 / 4.33				EN 455-2
Force at Break (N) During Aging	≥6 N				EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Concentration in µg/cm ²	Typical
Ammonium	0.055	Nitrate	0.277
Bromide	Not Detected	Nitrite	Not Detected
Calcium	0.813	Phosphate	Not Detected
Chloride	1.968	Potassium	0.095
Fluoride	0.005	Sodium	0.022
Lithium	Not Detected	Sulphate	0.538
Magnesium	0.020	Zinc	Not Detected

ORDERING INFORMATION

SIZE	S	M	L	XL
REORDER NO.	BLA301-S	BLA302-M	BLA303-L	BLA304-XL

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

