



Sterile Dispenser-Ready MCE Membrane Filter

Specifications

Membrane	Mixed Cellulose Ester (MCE)
Color	Black membrane with grid
Grid interval	3*3mm
Pore size	0.45μm
Diameter	47mm
Thickness	130±50μm
Water flow @bar*cm ² *min	≥70mL
Bubble point(kPa)	≥180
Bacterial retention capacity	LRV>7
Bacterial recovery rate	≥80%
pH range	pH 4-8
Max. Temp	130°C
Delivery condition	Sterile by Gamma Irradiation
Package	Dispenser Ready: 150pcs/box

Features

- High flow rate and flux
- Sterile individually packed and ready to use
- Grids for easy counting without restraining colony growth
- Excellent bacterial recovery rate
- Dispenser-ready membrane filter is compatible with other lab dispensers

Applications

- Microbial enumeration
- Bioburden testing of pharmaceutical and cosmetic products
- Bacterial analysis of water and beverages
- White membrane for colored bacteria while black membrane for light, white or transparent bacteria