

Aquanex De-Ionization packs includes right (A) and left (B) DI packs

The Thermo Scientific™ Aquanex™ is designed to truly help make life easier in lab with an intuitive touch screen interface, easy-access consumables, and smart dispensing features. The easy-to-use interface caters to users that need secured access and seamless event logging, system monitoring, and reporting. A new innovative cartridge design provides low TOC, a longer cartridge life, and an easy-to-change design.

Easily access consumables

- The Aquanex user interface provides step-by-step instructions on how to access and change consumables
- UV lamp, RO assembly, and ultrafilter consumables are easily accessed by removing the magnetic front panel
- No screws or complicated tools are required to change consumables

