

Herasafe 2030i 0.9m inc UV-C

The Herasafe 2030i series of biological safety cabinets leverages the legacy of BSC's global design leadership with innovations that optimize pollution prevention, user comfort and information management to deliver the most technologically advanced Class II cabinets available. Equipped with intelligent self-monitoring features, an intuitive touchscreen graphical interface, and built-in connectivity to preserve essential performance data, the HeraSafe 2030i system delivers unmatched precision performance for complex cell culture applications and other sensitive jobs that require product, personnel, and environmental safety to be ensured.

- **SmartFlow Plus:** This patented design performs automatic compensation that ensures perfect balance in airflows, regardless of the load condition of the protective HEPA filters and the clogging of the vents. This means that with a Thermo Scientific BSC, the level of containment is consistent every day, not just on the first day., Self-adjusting crossbeam UV light, Touchscreen GUI:
, All-316 stainless steel interior, LED Interior Lighting, Versatile SmartPort Access, Cloud-compatible data output, Optional ergonomic armrests

