

STER i160 CO2 SS TC 230V LK

Achieve your next breakthrough. The Thermo Scientific™ Forma™ Steri-Cycle™ i160 CO₂ 165 L Incubator is specifically designed to deliver the performance reliability, ease of operation, and value required to support a range of culturing needs from basic research to the demanding, leading-edge applications, so you're ready for whatever comes next. Better solutions for optimal cell growth. Our CO₂ incubators with non-corrosive stainless-steel or copper interior chambers provide the ideal in-vitro environment.

Direct Heat Chamber

- 165L electropolished stainless steel chamber with compact footprint, easily stackable
- Adjustable, perforated shelving
- Reversible exterior door for added flexibility

THRIVE Active Airflow

- In-chamber fan gently distributes clean, humidified air throughout the chamber ensuring homogeneous conditions
- Fast recovery of all parameters ten minutes or less for stable culturing conditions
- Auto-stop function disables fan operation during door openings to minimize air exchange

